

Purism Disposable Surgical Mask

Description:

Purism Disposable Surgical Mask is suitable for the daily protection of some non-oil particles, such as influenza virus, dust, pollen, etc., with good tightness, strong adsorption and good barrier.

Features:

Appearance:	plane-type
Size:	17.5cm×9.5cm

List of raw materials:

Polypropylene spunbonded non-woven fabric
Polypropylene melt-blown non-woven fabric
Polyethylene nose clip
Polyester / Polyurethane composites

Notes:

1. Please do not use the mask when the mask is contaminated or the whole mask is damaged and deformed.
2. This mask is disposable only and cannot be reused.
3. The side with the nose clip is facing outward when wearing.

Storage:

It is recommended to store in a clean, dry, ventilated and non-corrosive gas environment, away from ignition sources or combustible materials.

Important: the data in this paper is reliable and only for reference. We are not responsible for any results obtained by methods beyond our control. The user is responsible for determining the applicability of the product and taking appropriate measures to prevent possible loss and personal injury during storage and use of the product.